

This report was created by ACD/NMR Processor Academic Edition. For more information go to www.acdlabs.com/nmrproc/

Acquisition Time (sec)	3.1719	Comment	PROTON CDCl3 /opt/topspin GlassLab 32		Date	04 Jul 2013 10:27:28	
Date Stamp	04 Jul 2013 10:27:28						
File Name	C:\Users\Gregg Lowery\Documents\GlassLab\Summer 2013\NMR\Jul04-2013-GlassLab Pure p-toyl benzofulvene started 6.19.13\10\fid						
Frequency (MHz)	500.13	Nucleus	1H	Number of Transients	16	Origin	spect
Original Points Count	32768	Owner	nmr	Points Count	32768	Pulse Sequence	zg30
Receiver Gain	45.20	SW(cyclical) (Hz)	10330.58	Solvent	CHLOROFORM-d	Spectrum Offset (Hz)	3054.6028
Spectrum Type	STANDARD	Sweep Width (Hz)	10330.26	Temperature (degree C)	24.560		

